







The iGo2 rechargeable batteries can easily be exchanged without removing the unit from the carry case.









Accessories







iGo2 Charging Station (Battery not included)



DC Power Cord



Carrying Case



Rechargeable **Battery**

Call For Bundle Pricing

Battery Duration

Travel Bundles	Pulse Setting 2	Pulse Setting 3
+1 additional battery	7 Hours	6 Hours
+2 additional batteries	10.5 Hours	9 Hours







Comparisons Table

Manufacturer	Drive Devilbiss	Philips Respironics	GCE	Chart/Caire Airsep	Kingon
Product name and model	iGo2	SimplyGo Mini	Zen-O-lite	FreeStyle Comfort	P2
Device setting	1to 5	1to 5	1to 5	1to 5	1to 5
Max O2 Production (in ML)	1014	1000	1050	1050	1000
O ² Concentration	90% + 4/-3	87% - 96%	87% - 96%	90% + 5.5/-3	90% + 6/-3
Weight with Battery	2.2 kg	2.2 kg	2.5 kg	2.3-2.7 kg	1.97 kg
Sound Level at Pulse Setting 2	37dBA	42 dBA	37dBA	43 dBA	49 dBA 🜟
Conserver trigger sensitivity	0.05 cm H ² O	0.20 cm H²O	0.12 cm H ² O	0.5 cm H ² O	0.12 cm H²O
Battery Duration Pulse at Pulse Setting 2	Min 3.5 Hrs	Up to 4.5 Hrs	Up to 4 Hrs	Up to 4 Hrs	3 Hrs 20 minutes



A generator running at **50 dBA** is Twice as Noisy as a generator running at **40 dBA**.

Features



Revolutionary **SmartDose**TM

Technology automatically adjusts oxygen production & delivery according to patient activity.



Unmatched trigger sensitivity, meaning that the iGo2 will continue to work whilst you are sleeping.



Durable, **high-quality overmold design**provides an extra layer of protection against impact.



Federal Aviation Administration (FAA)
compliant for in-flight use.



The iGo2 POC can be operated whilst in its carry case. The over shoulder bag has been designed so that the device LCD screen and controls are easily accessible.



Covered by **medical aid** reimbursement.